	Gebrauchsanweisung	2
	Instructions for use	5
	Manuel d'utilisation.....	8
	Istruzioni per l'uso	11
	Instrucciones de uso.....	14
	Instrukcje użytkowania	17
	Instruktioner för användning	20
	Instruktioner til brug.....	23
	Bruksanvisning.....	26
	Käyttöohjeet	29
	Gebruiksaanwijzing	32
	Οδηγίες χρήσης.....	35
	Upute za uporabu.....	38
	Instrucțiuni de utilizare	41
	Naudojimo instrukcijos	44
	Návod k použití.....	47
	Návod na použitie	50
	Használati utasítás.....	53

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDKTLL

CopraTemp

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

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Produkt/ Produkttyp

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Instructions for use

Product/ Product type

PMMA blanks for the fabrication of temporary and permanent dentures and splints.

Material type

Polymethyl methacrylate (PMMA) - Class IIa medical device

Polymer-based crown and veneering material type 2 class 1 according to DIN EN ISO 10477 (all product variants)

Denture base acrylic type 1 class 1 according to DIN EN ISO 20795-1 (CopraTemp Base)

As the CAD/CAM blanks are polymerised during the manufacturing process, the polymerisation process at the user's premises is no longer necessary.

Users

Professional users who fabricate customised dental restorations/splints.

Target group

Patients of all ages and genders who require customised dental restorations.

Indication / Intended use

The products are medical devices intended for the manufacture of dental prostheses for long-term use, which are partially introduced into the human body through a clinical procedure, but are inserted into the teeth or attached to the tooth structure with suitable fixing material and are therefore classified as Class IIa.

CopraTemp and CopraTemp ML (coloured)

The tooth-coloured products are suitable for the fabrication of temporary crowns, bridges, inlays, onlays, etc., which are later replaced by the final prosthesis. However, they can also remain in the oral cavity as a permanent restoration.

CopraTemp (Clear)

The product is colourless and is used to make splints. Permanent retention in the oral cavity is indicated.

CopraTemp und CopraTemp ML (Base)

The pink-coloured product is suitable for the fabrication of denture bases, which are later replaced by the final denture. However, permanent retention in the oral cavity as a permanent restoration is also possible.

Contraindication

Do not use in case of proven hypersensitivity to one or more ingredients.

Do not use if there is insufficient space.

Veneering ceramics/ glaze

All ZrO₂ veneering ceramics/ glazes

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDKTLL

CopraTemp

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Material properties/ Technical data

Material composition: CopraTemp and CopraTemp ML coloured CopraTemp and CopraTemp ML Base	Polymethyl methacrylate, colour pigments
Material composition: CopraTemp Clear	Polymethyl methacrylate
Flexural strength (Weibull)	> 100 MPa
Water absorption	~ 20 to 23 µg/mm ³
Solubility	< 1 µg/mm ³
Vickers hardness	> 20 HV
E-modulus	> 2500 MPa
Residual monomer	< 1%

Processing

The desired restoration is milled from the selected blank. Acrylic milling cutters are recommended for processing. Polishing paste should only be used with soft polishing brushes. Please ensure that the material is not overheated.

A minimum wall thickness of 1.0 mm must be observed in the occlusal region and 0.6 mm in the cervical region. The connection between crowns or pontics must be at least 10 mm² in the anterior region and at least 15 mm² in the posterior region. No more than 2 pontics are permitted in the posterior region.

Veneering

All commercially available veneering materials that are suitable for PMMA - The respective instructions for use of the manufacturer of the bonder, veneering or luting material used must be observed.

Luting of plastic teeth

Plastic teeth can be attached to denture bases made of CopraTemp Base using a suitable luting material. The respective instructions for use of the manufacturer of the bonder, veneering or luting material used must be observed.

Safety instructions

Warning: Always wear respiratory protection (filter class FFP2), tight-fitting safety goggles, protective gloves and protective clothing and switch on extraction equipment. The dust generated during processing of this product can cause irritation to the skin/eyes/respiratory tract. Avoid inhalation/contact with skin, mouth and eyes. Do not eat or drink while working. Keep away from food and drink. Wash hands after handling. Remove contaminated clothing and protective equipment before entering eating areas. Keep away from sources of ignition. Do not smoke.

Storage

Store in a dry place. Protect from moisture. Store in the original packaging.

Waste disposal

Dispose of the product and packaging in accordance with local/ regional/ national/ international regulations.

Note

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the competent authority of the Member State in which the user and/or patient is established.

The Summary of Safety and Clinical Performance according to (EU) MDR 2017 / 745 Article 32 (SSCP) is available in the European database for medical devices (Eudamed) under the corresponding base UDI-DI as soon as the corresponding database module exists.

<https://ec.europa.eu/tools/eudamed>

SRN: DE-MF-000021120

Basis UDI-DI: ++EWPDCkTLL

CopraTemp

Whitepeaks Dental Solutions GmbH • Weikenrott 17 • 46499 Hamminkeln • Germany • +49 281 2064580

Explanation of the labelling on the packaging

Confirmation: The product complies with the applicable European directives

Catalogue number/ item number

Batch designation

Unique product identification

Number of products in the packaging

Please observe the instructions for use

www.wp-ifu.de/PMMA

Path to the electronic instructions for use

Instructions for use in paper form available within 7 days

Medical device

Date of manufacture

Manufacturer

Swiss authorised representative

Rx only

Caution: Federal law restricts this product to sale by or on the order of a licensed physician or dentist.

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